



## Trigno™ Avanti

Multi-purpose human movement sensor for the lab, classroom, and clinic

Trigno™ Avanti, the next generation wireless EMG and movement sensor system from Delsys, packs more technology than ever before while being our smallest sensor to date.

A true engineering feat, Trigno™ Avanti combines our patented EMG design with a high resolution inertial measurement unit (IMU) and flexible bandwidth. Designed for indoor and outdoor settings, Trigno™ Avanti includes an all new direct to tablet experience with research-oriented apps and tools.

Order Today



- EMG + IMU Sensor
- Direct to PC or Tablets
- Button-less Operation
- LED Status Indicators
- SDK/API Support for External Platforms
- Flexible EMG Bandwidth
- On-board Signal Processing & Data Storage
- Compatible with Existing Trigno™ Systems



[View this email in your browser](#)

You are receiving this email because of your relationship with Delsys. Please [reconfirm](#) your interest in receiving emails from us. If you do not wish to receive any more emails, you can [unsubscribe here](#).



Delsys Incorporated | 23 Strathmore Road | Natick, Massachusetts 01760

[Unsubscribe](#) | [Manage Subscription](#) | [Forward Email](#)



★ **This is a Test Email only.**

This message was sent for the sole purpose of testing a draft message.